


InstaMatch™

Version 1.0

Managing platelet refractory patients can be troublesome, time-consuming and costly. Finding appropriate blood products is essential to reduce a patient's risk of developing alloantibodies and to treat those who have become refractory.



Match every platelet every time

The InstaMatch web-based product matching and reservation software helps clinicians select the best-matched platelet product for a patient, based upon a grade and site-configurable compatibility score, and reserve it without logging into SafeTrace Tx®. An interface to a donor or transfusion management system¹ makes patient, donor and product information electronically available, automating the search and match process.

Compatible Products													
Unit ID	Product Code	Product Description	Grade	Score	Expiration Date/Time	Linked Patient MRN	Linked Patient Name	Antigens	ABO/Rh	Volume	CMV	Yield	Reserve/Cancel
W123416111506	E3045V00	APHERESIS PLATELETS ACD-A XX Z0-24C Irradiated	A	0	11/25/2016 11:59:00 PM	VAL002	Avalid, Jane	A1, A2, B7, B8	OP	350	Not Available	3.1x10 ¹¹	Reserve/Cancel
W123416111510	E3045V00	APHERESIS PLATELETS ACD-A XX Z0-24C Irradiated	B1U	0	11/25/2016 11:59:00 PM			A1, A2, B7	OP	150	Not Available	3.5x10 ¹¹	Reserve
W123416111507	E3045V00	APHERESIS PLATELETS ACD-A XX Z0-24C Irradiated	B2U	0	11/25/2016 11:59:00 PM	VAL002, VAL003	Avalid, Jane Bvalid, Jacob	A1, B7	OP	250	Not Available	1.5x10 ¹¹	Reserve/Cancel
W123416111509	E3045V00	APHERESIS PLATELETS ACD-A XX Z0-24C Irradiated	B1X	6	11/25/2016 11:59:00 PM			A1, A11, B7, B8	OP	250	Not Available	1.2x10 ¹¹	Reserve
W123416111508	E3045V00	APHERESIS PLATELETS ACD-A XX Z0-24C Irradiated	C	18	11/25/2016 11:59:00 PM			A1, A2, B27	OP	250	Not Available	2.1x10 ¹¹	Reserve

InstaMatch calculates a grade and site-configurable degree-of-compatibility score using a patient's HLA-type, blood attributes and special needs – so clinicians can select and reserve the best-matched product every time.

¹ Blood Establishment Computer Systems (BECS)

Using the InstaMatch™ product matching and reservation software you can:

Automate searching and matching

- Electronically receive patient and product information from a donor or transfusion management system
- Search for products based upon antigen avoidance and blood attributes (ABO-Rh, CVM status)
- Enable a consistent approach patient-to-patient, clinician-to-clinician

Manage refractory patients and help reduce alloantibodies risk

- Calculate a grade and site-configurable degree-of-compatibility score
- Select and reserve the best-matched platelet product

Streamline record keeping and coordinate user access

- Maintain search, edit and configuration change audit history
- Limit and control access with different user roles

Reservation Confirmation

You are about to reserve Unit ID: **W123416111506**
For Patient MRN: **VAL002**
Current Product Status: **Product in available inventory**

Expected Transfusion Date:
10/07/2016 1:29 PM
(UTC-05:00) Eastern Time (US & Canada) ▼

Comment
Canned Comment ▼
Patient arriving at 8:30 am]

Auto-reserve the best-matched product

Authorized clinicians can reserve a platelet product for a specific patient or cancel an existing reservation. Once a reservation is made, the product is electronically linked to the appropriate patient in SafeTrace Tx® - automating the process.

InstaMatch: Reserve/Cancel Unit

Unit Details

Unit ID: W123416111506
Product Code: E3045V00
Product Status: Product in available inventory
Description: APHERESIS PLATELETS|ACD-A|XX/20-24C|Irradiated
Location: ABCD
Expiration Date: 11/25/2016 11:59:00 PM
CMV: Not Available
Antigens: A1, A2, B7, B8
ABO/Rh: O P
Yield: 3.1x10¹¹
Comments: No comments

Search for a Patient

Q VAL002

✓ ABO/Rh Match: Compatible

	Grade	Score
Avalid, Jane VAL002	A	0

Service Provider: P48116
Date of Birth: 01/24/1966
Gender: Female
ABO/Rh: A P
Special Needs: None

✓ Reservation message sent. ✕

For more information or to schedule a demo, call Customer Service at 877.996.7877, email HSScustomerservice@haemonetics.com, or contact your Haemonetics Account Director.

For a list of worldwide office locations and contact information, visit www.haemonetics.com/officelocations

© 2016 Haemonetics Corporation. Haemonetics, InstaMatch, and SafeTrace Tx are trademarks or registered trademarks of Haemonetics Corporation in the USA, other countries, or both. 10.2016 COL-COPY-001077-US(AA)